

## Experiment No. 06


### Determination of Chemical Oxygen Demand (COD)

#### Objective:

The objective of this experiment is to determine the chemical oxygen demand (COD) of given water sample.

#### Chemical Oxygen Demand:

The amount of oxygen required to oxidize all the organic matter present in wastewater is known as chemical oxygen demand (COD) and is expressed in terms of milligrams of oxygen required per litre of water, mg/l. The COD test allows measurement of oxygen demand of the waste in terms of the total quantity of oxygen required for oxidation of the waste to carbon dioxide and water. The test is based on the fact that all organic compounds, with a few exceptions, can be oxidized chemically by the action of strong oxidizing agents under acidic condition.



The reaction in equation (1) involves conversion of organic matter to carbon dioxide and water regardless of the biological assimilability of the substance. For example, glucose and lignin (biologically inert substance) are both oxidized completely by the chemical oxidant. As a result, COD values are greater than BOD values, especially when biologically resistant organic matter is present.

The main advantage of COD test is the short time required for evaluation. The determination can be made in about 3 hours rather than the 5-days required for the measurement of BOD. For this reason, it is used as a substitute for the BOD test in many instances.

One of the major limitations of COD test is its inability to differentiate between biodegradable and non-biodegradable organic matter. In addition, it does not provide any evidence of the rate at which the biologically active material would be stabilized under conditions that exist in nature.

#### Environmental significance:

COD is often measured as a rapid indicator of organic pollutant in water. It is normally measured in both municipal and industrial wastewater treatment plants and gives an indication of the efficiency of the treatment process.

*According to the Environment Conservation Rules (1997), drinking water standard for COD is 4.0 mg/l. For wastewater effluent allowable concentration of COD varies from 200-400*

*mg/l depending on discharge point of the effluent (e.g., inland water, irrigation land, public sewer etc.).*

### **Theory on experimental method:**

Potassium dichromate or potassium permanganate is usually used as the oxidizing agent in the determination of COD. In this test, potassium permanganate would be used. Potassium permanganate is selective in the reaction and attacks the carbonaceous and not the nitrogenous matter.

In any method of measuring COD, an excess of oxidizing agent must be present to ensure that all organic matter is oxidized as completely as possible within the power of the reagent. This requires that a reasonable excess be present in all samples. It is necessary, therefore, to measure the excess in some manner so that the actual amount can be determined. For this purpose, a solution of a reducing agent (e.g., ammonium oxalate) is usually used.

### **Reagents:**

1. Diluted sulfuric acid
2. Standard potassium permanganate solution
3. Standard ammonium oxalate solution

### **Apparatus:**

1. Beaker (250 ml)
2. Dropper
3. Stirrer

### **Procedure:**

1. Pipette 100 ml of the sample into a 250 ml Erlenmeyer flask.
2. Add 10 ml of diluted sulfuric acid and 10 ml of standard  $\text{KMnO}_4$  solution.
3. Heat the flask in a boiling water bath for exactly 30 minutes, keeping the water in the bath above the level of the solution in the flask. The heating enhances the rate of oxidation reaction in the flask.
4. If the solution becomes faintly coloured, it means that most of the potassium permanganate has been utilized in the oxidation of organic matter. In such case, repeat the above using a smaller sample diluted to 100 ml with distilled water.
5. After 30 minutes in the water bath, add 10 ml of standard ammonium oxalate solution into the flask. This 10 ml ammonium oxalate, a reducing agent, is just equivalent to the 10 ml potassium permanganate (oxidizing agent) added earlier. The excess of reducing agent now remaining in the flask is just equivalent to the amount of the oxidizing agent used in the oxidation of organic matter.

6. The quantity of ammonium oxalate remaining in the flask is now determined by titration with standard potassium permanganate. Titrate the content of the flask while hot with standard potassium permanganate to the first pink coloration. Record the ml of potassium permanganate used.

**Observation and Calculation:**

$$\text{COD (mg/l)} = (\text{ml of KMnO}_4 \text{ used in step 6} \times 100) / \text{ml of sample used}$$

**Result:**

**Comments on Result:**